



February 1, 2022

Dear Vesco Medical Customers & Colleagues,

In 2020, The Global Enteral Device Supplier Association (GEDSA) released an ENFit Connector Conversion Schedule in which members agreed to no longer manufacture legacy enteral feeding devices and transition connectors. Due to the persistence of the pandemic and related global supply chain challenges, implementation timelines for most manufacturers have been delayed. However, our resolve to maximize patient safety by reducing the risk of enteral tubing misconnections remains. Now is the time to transition to safer ISO 80369-3 standard ENFit connectors once and for all.

To support a cohesive and safe transition, Vesco will continue to offer products with ENFit Transition Connectors (TCs) to our committed customers for a period of time after the last anticipated manufactured date while supplies last. See adjacent table with crosswalk of existing items with TCs and the corresponding item numbers without TCs. Separate items numbers should improve tracking inventory levels and ordering with or without TCs. Individually packaged Sterile, Clean Non-Sterile and bulk, non-sterile **ENFit syringes without TCs are readily available today**. The following outlines our projected availability and the month we anticipate running out of product with TCs in pack:

Product Line	Last Manufactured Date
Gravity Sets with Transition Connectors	July, 2022
ENFit Syringes with Transition Connectors	July, 2022
ENFit Transition Connectors (sold separately)	September, 2022

At Vesco we are **All ENFit All the Time!** We remain fully committed to supporting the industry-wide transition. We carry a full range of ENFit products to support your evolving enteral device needs. Patient safety is of the utmost importance. Our team of ENFit experts are prepared to help solve any transition challenges or questions. For assistance, please reach out to us at (855) 239-0131 or info@vescomedical.com.

Thank you for your attention to this important matter.

Kind Regards,

Tom Hancock
President & Founder, Vesco Medical

¹ To read GEDSA's full communication and/or learn more about the industry's recommendations to support a smooth ENFit transition, please visit www.stayconnected.org or contact info@gedsa.org.

Vesco Medical



Items with/without ENFit Transition Connectors

Item Description	Transition Set (w/TC) Item #	Final Set (no TC) Item #	Packaging
40mm Screw Cap Gravity Set with ENFit	VED-020	VED-022	30/Case
ENPlus Spike Gravity Set with ENFit	VED-026	VED-028	30/Case
500mL Gravity Bag Set (LL) with ENFit	VED-038	VED-040	30/Case
1000mL Gravity Bag Set (LL) with ENFit	VED-032	VED-034	30/Case
1000mL Gravity TOP FILL Bag Set (PVC) with ENFit	VED-045	VED-046	30/Case
1000mL Gravity LARGE BORE Bag Set with ENFit	VED-048	VED-049	30/Case
5mL ENFit Tip Syringe with TC	VED-605TC	VED-605	500/Case
10mL ENFit Tip Syringe with TC	VED-610TC	VED-610	500/Case
20mL ENFit Tip Syringe with TC	VED-620TC	VED-620	500/Case
35mL ENFit Tip Syringe with TC	VED-635TC	VED-635	500/Case
60mL ENFit Tip Syringe with TC	VED-660TC	VED-660	180/Case
ENFit Transition Connector	VED-001	N/A	10/Bag x 30 Bags = 300/Case